



AHCY rabbit pAb

货号	YP-Ab-18347
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
靶点	
基因名称	AHCY SAHH
蛋白名称	Adenosylhomocysteinase (AdoHcyase) (S-adenosyl-L-homocysteine hydrolase)
免疫原	Synthesized peptide derived from human AHCY
特异性	This antibody detects endogenous levels of AHCY at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	48kD
功能	Catalyzes the hydrolysis of S-adenosyl-L-homocysteine to form adenosine and homocysteine . Binds copper ions (By similarity).
细胞定位	Cytoplasm . Melanosome . Identified by mass spectrometry in melanosome fractions from stage I to stage IV.
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images