



CHCHD6 rabbit pAb

货号	YP-Ab-18353
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
靶点	
基因名称	CHCHD6 CHCM1
蛋白名称	Coiled-coil-helix-coiled-coil-helix domain-containing protein 6, mitochondrial (Coiled coil helix cristae morphology protein 1)
免疫原	
特异性	This antibody detects endogenous levels of CHCHD6 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	26kD
功能	Component of the MICOS complex, a large protein complex of the mitochondrial inner membrane that plays crucial roles in the maintenance of crista junctions, inner membrane architecture, and formation of contact sites to the outer membrane.
细胞定位	Mitochondrion inner membrane ; Lipid-anchor . Mitochondrion .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

