



ANXA9 mouse mAb

货号	YP-mAb-18318
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
基因名称	ANXA9 ANX31
蛋白名称	Annexin A9 (Annexin XXXI) (Annexin-31) (Annexin-9) (Pemphaxin)
免疫原	Synthesized peptide derived from human ANXA9
特异性	This antibody detects endogenous levels of ANXA9 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	38kD
功能	Low affinity receptor for acetylcholine known to be targeted by disease-causing pemphigus vulgaris antibodies in keratinocytes.
组织表达	Expressed in the stratified squamous skin epithelium, but not in epithelia of other types (at protein level).
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images