



## RAE1 mouse mAb

货号	YP-mAb-18322
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
靶点	
基因名称	RAE1 MRNP41
蛋白名称	mRNA export factor (Rae1 protein homolog) (mRNA-associated protein mrnp 41)
免疫原	Synthesized peptide derived from human RAE1
特异性	This antibody detects endogenous levels of RAE1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	40kD
功能	Acts as mRNA export factor involved in nucleocytoplasmic transport . Plays a role in mitotic bipolar spindle formation . May function in attaching cytoplasmic mRNPs to the cytoskeleton both directly or indirectly .
细胞定位	Cytoplasm . Nucleus . Cytoplasm, cytoskeleton, spindle pole . Nucleus envelope . Recruited from interphase nuclei to spindle MTs during mitosis. .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

