



LRRTM3 mouse mAb

货号	YP-mAb-18326
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
基因名称	LRRTM3 UNQ803/PRO1693
蛋白名称	Leucine-rich repeat transmembrane neuronal protein 3
免疫原	Synthesized peptide derived from human LRRTM3
特异性	This antibody detects endogenous levels of LRRTM3 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	64kD
功能	Exhibits a limited synaptogenic activity in vitro, restricted to excitatory presynaptic differentiation (By similarity). May play a role in the development and maintenance of the vertebrate nervous system.
细胞定位	Cell membrane ; Single-pass type I membrane protein . Cell junction, synapse, postsynaptic cell membrane ; Single-pass type I membrane protein .
组织表达	Expressed in neuronal tissues.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images