



DCTN6 mouse mAb

货号	YP-mAb-18331
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
基因名称	DCTN6 WS3
蛋白名称	Dynactin subunit 6 (Dynactin subunit p27) (Protein WS-3)
免疫原	Synthesized peptide derived from human DCTN6
特异性	This antibody detects endogenous levels of DCTN6 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	21kD
功能	
细胞定位	Cytoplasm, cytoskeleton . Chromosome, centromere, kinetochore .
组织表达	Ubiquitous.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images