



## P4HA2 mouse mAb

货号	YP-mAb-18339
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
基因名称	P4HA2 UNQ290/PRO330
蛋白名称	Prolyl 4-hydroxylase subunit alpha-2 (4-PH alpha-2) (Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2)
免疫原	Synthesized peptide derived from human P4HA2
特异性	This antibody detects endogenous levels of P4HA2 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	59kD
功能	Catalyzes the post-translational formation of 4-hydroxyproline in -Xaa-Pro-Gly-sequences in collagens and other proteins.
细胞定位	Endoplasmic reticulum lumen.
组织表达	Expressed in the heart, placenta, lung and pancreas.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

### Products Images