



CALCOCO1 rabbit pAb

货号	YP-Ab-18354
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
靶点	
基因名称	CALCOCO1 KIAA1536 PP13275 UNQ2436/PRO4996
蛋白名称	Calcium-binding and coiled-coil domain-containing protein 1 (Calphoglin) (Coiled-coil coactivator protein) (Sarcoma antigen NY-SAR-3)
免疫原	Synthesized peptide derived from human CALCOCO1
特异性	This antibody detects endogenous levels of CALCOCO1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	76kD
功能	Functions as a coactivator for aryl hydrocarbon and nuclear receptors (NR). Recruited to promoters through its contact with the N-terminal basic helix-loop-helix-Per-Arnt-Sim (PAS) domain of transcription factors or coactivators, such as NCOA2. During ER-activation acts synergistically in combination with other NCOA2-binding proteins, such as EP300, CREBBP and CARM1. Involved in the transcriptional activation of target genes in the Wnt/CTNNB1 pathway. Functions as a secondary coactivator in LEF1-mediated transcriptional activation via its interaction with CTNNB1. Coactivator function for nuclear receptors and LEF1/CTNNB1 involves differential utilization of two different activation regions (By similarity). In association with CCAR1 enhances GATA1- and MED1-mediated transcriptional activation from the gamma-globin promoter during erythroid differentiation of K562 erythroleukemia cells . ; Seems to enhance inorganic pyrophosphatase thus activating phosphoglucomutase (PMG). Probably functions as component of the calphoglin complex, which is involved in linking cellular metabolism (phosphate and glucose metabolism) with other core functions including protein synthesis and degradation, calcium signaling and cell growth.
细胞定位	Cytoplasm. Nucleus. Shuttles between nucleus and cytoplasm. .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)



有关注意事项

Avoid repeated freezing and thawing!

使用建议

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images