



## LRRC17 rabbit pAb

货号	YP-Ab-18356
同位型	IgG
应用	WB
种属	Human, Mouse,Rat
靶点	
基因名称	LRRC17 P37NB UNQ3076/PRO9909
蛋白名称	Leucine-rich repeat-containing protein 17 (p37NB)
免疫原	Synthesized peptide derived from human LRRC17
特异性	This antibody detects endogenous levels of LRRC17 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	49kD
功能	Involved in bone homeostasis. Acts as a negative regulator of RANKL-induced osteoclast precursor differentiation from bone marrow precursors (By similarity).
细胞定位	Secreted, extracellular space .
组织表达	Expressed in osteoblast cell lines. Well expressed in ovary, heart, pancreas, skeletal muscle, lung, and fetal kidney and lung and only at the basal levels in the other tissues examined including adult kidney. More expressed in S-type neuroblastoma cells than in N-type neuroblastoma cells.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

