



DTX2 rabbit pAb

货号	YP-Ab-18377
同位型	IgG
应用	WB
种属	Human,Mouse
基因名称	DTX2 KIAA1528 RNF58
蛋白名称	Probable E3 ubiquitin-protein ligase DTX2 (Protein deltex-2) (Deltex2) (hDTX2) (RING finger protein 58)
免疫原	Synthesized peptide derived from human DTX2
特异性	This antibody detects endogenous levels of DTX2 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	68kD
功能	Regulator of Notch signaling, a signaling pathway involved in cell-cell communications that regulates a broad spectrum of cell-fate determinations. Probably acts both as a positive and negative regulator of Notch, depending on the developmental and cell context. Mediates the antineural activity of Notch, possibly by inhibiting the transcriptional activation mediated by MATCH1. Functions as a ubiquitin ligase protein in vitro, suggesting that it may regulate the Notch pathway via some ubiquitin ligase activity.
细胞定位	Cytoplasm . Nucleus . Predominantly cytoplasmic. Partially nuclear. .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images