



RAB3B/C rabbit pAb

Catalog No	YP-Ab-18409
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RAB3B
Protein Name	Ras-related protein Rab-3B
Immunogen	Synthesized peptide derived from human RAB3B/C
Specificity	This antibody detects endogenous levels of RAB3B/C at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	24kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side . Golgi apparatus . Colocalizes with GAS8/DRC4 in the Golgi apparatus. .
Tissue Specificity	
Function	Protein transport. Probably involved in vesicular traffic (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

