



## KLC2 Mouse mAb

<b>Catalog No</b>	YP-mAb-18471
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse,Rat
<b>Applications</b>	WB
<b>Gene Name</b>	KLC2
<b>Protein Name</b>	Kinesin light chain 2 (KLC 2)
<b>Immunogen</b>	Synthesized peptide derived from human KLC2
<b>Specificity</b>	This antibody detects endogenous levels of KLC2 at Human
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	68kD
<b>Cell Pathway</b>	Cytoplasm, cytoskeleton . Lysosome membrane ; Peripheral membrane protein ; Cytoplasmic side .
<b>Tissue Specificity</b>	
<b>Function</b>	Kinesin is a microtubule-associated force-producing protein that plays a role in organelle transport. The light chain functions in coupling of cargo to the heavy chain or in the modulation of its ATPase activity (Probable). Through binding with PLEKHM2 and ARL8B, recruits kinesin-1 to lysosomes and hence direct lysosomes movement toward microtubule plus ends .
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images