



# Zip9 Mouse mAb(Mix-mA)

<b>Catalog No</b>	YP-Ab-17216
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human, Mouse,Rat
<b>Applications</b>	WB,IHC
<b>Gene Name</b>	SLC39A9 ZIP9 UNQ714/PRO1377
<b>Protein Name</b>	Zip9
<b>Immunogen</b>	Synthesized peptide derived from human Zip9
<b>Specificity</b>	This antibody detects endogenous levels of Zip9 at Human, Mouse,Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Mouse,monoclonal
<b>Purification</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000 IHC 1:50-200
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Zinc transporter ZIP9 (Solute carrier family 39 member 9) (Zrt- and Irt-like protein 9) (ZIP-9)
<b>Observed Band</b>	33kD
<b>Cell Pathway</b>	Membrane ; Multi-pass membrane protein .
<b>Tissue Specificity</b>	
<b>Function</b>	function:May act as a zinc-influx transporter.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.,
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

