



Plexin A1 Rabbit mAb

Catalog No	YP-Ab-17768
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,ICC/IF,FC
Gene Name	PLXNA1
Alternative Names	NOV; NOV; Plexin-A1; Plexin1; PlexinA1; PLXN1; Plxna1
Research Field	Neuroscience
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 211 kDa; Observed MW: 211 kDa
Clonality	Monoclonal Antibody
Clonality No.	R08-2J4
Dilution	WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100
Immunogen	A synthesized peptide derived from human Plexin A1
Purification	Affinity Chromatography
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Plays a role in axon guidance, invasive growth and cell migration.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

