



TAX1BP1 Rabbit mAb

Catalog No	YP-Ab-17782
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB,ICC/IF
Gene Name	TAX1BP1
Alternative Names	CALCOCO3; D6Ertd404e; D6Ertd772e; PRO0105; T6BP; TAX1BP1; tax1bp1b; TXBP151
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 91 kDa; Observed MW: 91 kDa
Clonality	Monoclonal Antibody
Clonality No.	R07-1J3
Dilution	WB: 1/500-1/1000 IF: 1/50-1/200
Immunogen	A synthesized peptide derived from human TRAF6BP
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degraded by caspase-3-like family proteins upon TNF-induced apoptosis. May also play a role in the pro-inflammatory cytokine IL-1 signaling cascade.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

