



# ACAA1 rabbit pAb

货号	YP-Ab-18205
同位型	IgG
应用	WB
种属	Human
靶点	ACAA1
基因名称	ACAA1 ACAA PTHIO
蛋白名称	3-ketoacyl-CoA thiolase, peroxisomal (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase) (Peroxisomal 3-oxoacyl-CoA thiolase)
免疫原	Synthesized peptide derived from human ACAA1
特异性	This antibody detects endogenous levels of ACAA1 at Human
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	47kD
功能	Responsible for the thiolytic cleavage of straight chain 3-keto fatty acyl-CoAs (3-oxoacyl-CoAs) (Probable). Plays an important role in fatty acid peroxisomal beta-oxidation (Probable). Catalyzes the cleavage of short, medium, long, and very long straight chain 3-oxoacyl-CoAs (Probable).
细胞定位	Peroxisome .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images