



SNRPB2 rabbit pAb

货号	YP-Ab-18215
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	SNRPB2
基因名称	SNRPB2
蛋白名称	U2 small nuclear ribonucleoprotein B" (U2 snRNP B")
免疫原	Synthesized peptide derived from human SNRPB2
特异性	This antibody detects endogenous levels of SNRPB2 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	25kD
功能	Involved in pre-mRNA splicing as component of the spliceosome . Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA .
细胞定位	Nucleus .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images