



VAMP3/Cellubrevin rabbit pAb

货号	YP-Ab-18220
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	VAMP3
基因名称	VAMP3 SYB3
蛋白名称	Vesicle-associated membrane protein 3 (VAMP-3) (Cellubrevin) (CEB) (Synaptobrevin-3)
免疫原	Synthesized peptide derived from human VAMP3/Cellubrevin
特异性	This antibody detects endogenous levels of VAMP3/Cellubrevin at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	11kD
功能	SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.
细胞定位	Membrane ; Single-pass type IV membrane protein . Cell junction, synapse, synaptosome .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images