



# MIPEP rabbit pAb

货号	YP-Ab-18223
同位型	IgG
应用	WB
种属	Human;Mouse;Ra
靶点	MIPEP
基因名称	MIPEP MIP
蛋白名称	Mitochondrial intermediate peptidase (MIP) (EC 3.4.24.59)
免疫原	Synthesized peptide derived from human MIPEP
特异性	This antibody detects endogenous levels of MIPEP at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	78kD
功能	Cleaves proteins, imported into the mitochondrion, to their mature size.
细胞定位	Mitochondrion matrix.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images