



# CCPG1 rabbit pAb

货号	YP-Ab-18251
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	CCPG1
基因名称	CCPG1 CCP8 CPR8 KIAA1254
蛋白名称	Cell cycle progression protein 1 (Cell cycle progression restoration protein 8)
免疫原	Synthesized peptide derived from human CCPG1
特异性	This antibody detects endogenous levels of CCPG1 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	83kD
功能	Acts as an assembly platform for Rho protein signaling complexes. Limits guanine nucleotide exchange activity of MCF2L toward RHOA, which results in an inhibition of both its transcriptional activation ability and its transforming activity. Does not inhibit activity of MCF2L toward CDC42, or activity of MCF2 toward either RHOA or CDC42 (By similarity). May be involved in cell cycle regulation.
细胞定位	Cytoplasmic granule membrane ; Single-pass type II membrane protein .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

