



GOLGA3 rabbit pAb

货号	YP-Ab-18255
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	GOLGA3
基因名称	GOLGA3
蛋白名称	Golgin subfamily A member 3 (Golgi complex-associated protein of 170 kDa) (GCP170) (Golgin-160)
免疫原	Synthesized peptide derived from human GOLGA3
特异性	This antibody detects endogenous levels of GOLGA3 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	165kD
功能	Golgi auto-antigen; probably involved in maintaining Golgi structure.
细胞定位	Cytoplasm. Golgi apparatus, Golgi stack membrane; Peripheral membrane protein.
组织表达	Expressed in all tissues tested. Expressed in liver, testis, lung, heart, salivary gland and kidney.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

