

Sam50 rabbit pAb

货号	YP-Ab-18270
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	SAMM50
基因名称	SAMM50 CGI-51 SAM50 TRG3
蛋白名称	Sorting and assembly machinery component 50 homolog (Transformation-related gene 3 protein) (TRG-3)
免疫原	Synthesized peptide derived from human Sam50
特异性	This antibody detects endogenous levels of Sam50 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	52kD
功能	Plays a crucial role in the maintenance of the structure of mitochondrial cristae and the proper assembly of the mitochondrial respiratory chain complexes . Required for the assembly of TOMM40 into the TOM complex .
细胞定位	Mitochondrion outer membrane ; Multi-pass membrane protein. Cytoplasm . Mitochondrion .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

