



## NAT1 rabbit pAb

货号	YP-Ab-18280
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	NAT1
基因名称	NAT1 AAC1
蛋白名称	Arylamine N-acetyltransferase 1 (EC 2.3.1.5) (Arylamide acetylase 1) (Monomorphic arylamine N-acetyltransferase) (MNAT) (N-acetyltransferase type 1) (NAT-1)
免疫原	Synthesized peptide derived from human NAT1
特异性	This antibody detects endogenous levels of NAT1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	32kD
功能	Participates in the detoxification of a plethora of hydrazine and arylamine drugs. Catalyzes the N- or O-acetylation of various arylamine and heterocyclic amine substrates and is able to bioactivate several known carcinogens.
细胞定位	Cytoplasm.
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

