



TSEPA rabbit pAb

货号	YP-Ab-18284
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	THEMIS
基因名称	THEMIS C6orf190 C6orf207
蛋白名称	Protein THEMIS (Thymocyte-expressed molecule involved in selection)
免疫原	Synthesized peptide derived from human TSEPA
特异性	This antibody detects endogenous levels of TSEPA at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	71kD
功能	Plays a central role in late thymocyte development by controlling both positive and negative T-cell selection. Required to sustain and/or integrate signals required for proper lineage commitment and maturation of T-cells. Regulates T-cell development through T-cell antigen receptor (TCR) signaling and in particular through the regulation of calcium influx and phosphorylation of Erk.
细胞定位	Cytoplasm . Nucleus .
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

