



CKAP4 rabbit pAb

货号	YP-Ab-18288
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	CKAP4
基因名称	CKAP4
蛋白名称	Cytoskeleton-associated protein 4 (63-kDa cytoskeleton-linking membrane protein) (Climp-63) (p63)
免疫原	Synthesized peptide derived from human CKAP4
特异性	This antibody detects endogenous levels of CKAP4 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	66kD
功能	Mediates the anchoring of the endoplasmic reticulum to microtubules. ; High-affinity epithelial cell surface receptor for the FZD8-related low molecular weight sialoglycopeptide APF/antiproliferative factor. Mediates the APF antiproliferative signaling within cells.
细胞定位	Endoplasmic reticulum membrane ; Single-pass type II membrane protein. Cell membrane ; Single-pass type II membrane protein. Cytoplasm, cytoskeleton. Cytoplasm, perinuclear region. Translocates to the perinuclear region upon APF-stimulation.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

