



# ATOX1 rabbit pAb

货号	YP-Ab-18290
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	ATOX1
基因名称	ATOX1 HAH1
蛋白名称	Copper transport protein ATOX1 (Metal transport protein ATX1)
免疫原	Synthesized peptide derived from human ATOX1
特异性	This antibody detects endogenous levels of ATOX1 at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	7kD
功能	Binds and deliver cytosolic copper to the copper ATPase proteins. May be important in cellular antioxidant defense.
组织表达	Ubiquitous.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

## Products Images