

DDX41 rabbit pAb

货号	YP-Ab-18291
同位型	IgG
应用	WB
种属	Human;Mouse
靶点	DDX41
基因名称	DDX41 ABS
蛋白名称	Probable ATP-dependent RNA helicase DDX41 (EC 3.6.4.13) (DEAD box protein 41) (DEAD box protein abstract homolog)
免疫原	Synthesized peptide derived from human DDX41
特异性	This antibody detects endogenous levels of DDX41 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	68kD
功能	Probable ATP-dependent RNA helicase. Is required during post-transcriptional gene expression. May be involved in pre-mRNA splicing.
细胞定位	Nucleus .
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images