



PIK3CA H1047R (ABT453) Mouse mAb


Catalog No	YP-Ab-15691
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	PIK3CA
Protein Name	5-bisphosphate 3-kinase 110 kDa catalytic subunit alpha; 5-bisphosphate 3-kinase catalytic subunit alpha isoform; caPI3K; CLOVE; CWS5; MCAP; MCM; MCMTC; MGC142161; MGC142163; p110 alpha; p110alpha; Phosphatidylin
Immunogen	Synthesized peptide derived from human PIK3CA H1047R
Specificity	The antibody can specifically recognize the h1047r mutation of human PIK3CA protein, and does not react with wild-type PIK3CA protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG2b, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: 500-1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5-bisphosphate 3-kinase 110 kDa catalytic subunit alpha; 5-bisphosphate 3-kinase catalytic subunit alpha isoform; caPI3K; CLOVE; CWS5; MCAP; MCM; MCMTC; MGC142161; MGC142163; p110 alpha; p110alpha; Phosphatidylinositol 3 kinase catalytic alpha polypeptide; Phosphatidylinositol 3 kinase catalytic 110 KD alpha; Phosphatidylinositol 4 5 bisphosphate 3 kinase catalytic subunit alpha; Phosphatidylinositol 4 5 bisphosphate 3 kinase catalytic subunit alpha isoform; Phosphatidylinositol 4,5 bisphosphate 3 kinase 110 kDa catalytic subunit alpha; Phosphatidylinositol-4; Phosphoinositide 3 kinase catalytic alpha polypeptide; PI 3 Kinase catalytic subunit alpha; PI3 kinase p110 subunit alpha; PI3-kinase subunit alpha; PI3K; PI3K-alpha; PI3KC A; PIK3C A; Pik3ca; PK3CA; PK3CA_HUMAN; PtdIns 3 kinase p110; PtdIns-3-kinase subunit alpha; PtdIns-3-kinase subunit p110-alpha; Serine/threonine protein kinase PIK3CA
Observed Band	




Cell Pathway	Cytoplasmic
Tissue Specificity	Gastric adenocarcinoma with PIK3CA H1047R protein expression
Function	<p>catalytic activity:ATP + 1-phosphatidyl-1D-myo-inositol 4,5-bisphosphate = ADP + 1-phosphatidyl-1D-myo-inositol 3,4,5-trisphosphate, .disease:Defects in PIK3CA are associated with breast cancer [MIM:114480],.disease:Defects in PIK3CA are associated with colorectal cancer (CRC) [MIM:114500],.disease:Defects in PIK3CA are associated with ovarian cancer [MIM:167000]. Ovarian cancer is the leading cause of death from gynecologic malignancy. It is characterized by advanced presentation with loco-regional dissemination in the peritoneal cavity and the rare incidence of visceral metastases. These typical features relate to the biology of the disease, which is a principal determinant of outcome.,disease:Defects in PIK3CA may underlie hepatocellular carcinoma (HCC) [MIM:114550],.disease:PI3KCA mutations affecting exons 9 and 20 display gender-and tissue-specific patterns, thus suggesting that the</p>
Background	<p>Phosphatidylinositol 3-kinase is composed of an 85 kDa regulatory subunit and a 110 kDa catalytic subunit. The protein encoded by this gene represents the catalytic subunit, which uses ATP to phosphorylate PtdIns, PtdIns4P and PtdIns(4,5)P2. This gene has been found to be oncogenic and has been implicated in cervical cancers. A pseudogene of this gene has been defined on chromosome 22. [provided by RefSeq, Apr 2016],</p>
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Human gastric adenocarcinoma tissue was stained with Anti-PIK3CA H1047R (ABT453) Antibody




Human gastric adenocarcinoma tissue was stained with Anti-PIK3CA H1047R (ABT453) Antibody



Human gastric adenocarcinoma tissue was stained with Anti-PIK3CA H1047R (ABT453) Antibody



Human gastric adenocarcinoma tissue was stained with Anti-PIK3CA H1047R (ABT453) Antibody



Human gastric adenocarcinoma tissue was stained with Anti-PIK3CA H1047R (ABT453) Antibody