



Caspase-10 B/C Polyclonal Antibody

Catalog No	YP-Ab-00342
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC-p;IF/ICC;ELISA
Gene Name	CASP10
Protein Name	Caspase10
Immunogen	The antiserum was produced against synthesized peptide derived from human Caspase 10. AA range:430-479
Specificity	Caspase-10 B/C Polyclonal Antibody detects endogenous levels of Caspase-10 B/C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALPS2; apoptotic protease Mch-4; CASP-10; CASPA; caspase-10; FADD-like ICE2; FAS-associated death domain protein interleukin-1B-convert
Observed Band	60
Cell Pathway	cytosol,CD95 death-inducing signaling complex,riposome,
Tissue Specificity	Detectable in most tissues. Lowest expression is seen in brain, kidney, prostate, testis and colon.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images