



IRS-1 (phospho-Ser612) rabbit pAb

Catalog No	YP-Ab-03629
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IF/ICC;ELISA;IHC-p
Gene Name	IRS1
Protein Name	IRS-1 (Ser612)
Immunogen	Synthesized phospho peptide around human IRS-1 (Ser612)
Specificity	This antibody detects endogenous levels of Human Mouse Rat IRS-1 (phospho-Ser612)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IF/ICC 1:50-200;ELISA 1:2000-20000;IHC-p 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Insulin receptor substrate 1 (IRS-1)
Observed Band	170
Cell Pathway	nucleus,cytoplasm,cytosol,plasma membrane,insulin receptor complex,caveola,intracellular membrane-bounded organelle,
Tissue Specificity	Epithelium, Eye, Skeletal muscle,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

