



Myc-Tag Monoclonal Antibody(5D6)

Catalog No	YP-Ab-04743
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	Myc Tag
Protein Name	Myc-Tag
Immunogen	Synthetic peptide EQKLISEEDL coupled to KLH.
Specificity	Antibody can detect C-Myc fusion proteins.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000 -10000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

