



GFP/eGFP-Tag mouse mAb

Catalog No	YP-Ab-04514
Isotype	IgG
Reactivity	EGFP;EYFP
Applications	WB;IF/ICC;
Gene Name	egfp
Protein Name	GFP Tag
Immunogen	Purified recombinant EYFP.tag full length expressed in E.coli.
Specificity	Transfected Only.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb dilution 1:1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	According to the fused protein
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images