



# ECFP-Tag Monoclonal Antibody(8F6), HRP Conjugated

|                                  |   |
|----------------------------------|---|
| <b>Catalog No</b>                | YP-Ab-04557   |
| <b>Isotype</b>                   | IgG   |
| <b>Reactivity</b>                | Species independent   |
| <b>Applications</b>              | WB  |
| <b>Gene Name</b>                 |   |
| <b>Protein Name</b>              |   |
| <b>Immunogen</b>                 |   |
| <b>Specificity</b>               | ECFP-Tag Monoclonal Antibody(8F6) HRP Conjugated, specially designed for your Western blot analysis.  |
| <b>Formulation</b>               | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.  |
| <b>Source</b>                    | Monoclonal, Mouse IgG   |
| <b>Purification</b>              | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.                              |
| <b>Dilution</b>                  | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:2000-5000. |
| <b>Concentration</b>             | 1 mg/ml   |
| <b>Purity</b>                    | ≥90%  |
| <b>Storage Stability</b>         | -20°C/1 year  |
| <b>Synonyms</b>                  | ECFP  |
| <b>Observed Band</b>             |   |
| <b>Cell Pathway</b>              |   |
| <b>Tissue Specificity</b>        |   |
| <b>Function</b>                  |   |
| <b>Background</b>                |   |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!  |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.               |



## Products Images