




SRT-Tag Monoclonal Antibody(4H9)

Catalog No	YP-Ab-04757
Isotype	IgG
Reactivity	Species independent
Applications	WB
Gene Name	
Protein Name	
Immunogen	Recombinant Protein of SRT-Tag
Specificity	The antibody detects SRT Tag recominant protein.
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB: 1:5000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	The SRT-tag is a peptide epitope comprised of 10 amino acids derived from the crystalline surface layer protein of Rickettsia typhi.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of SRT-Recombinant protein,
diluted at 1) 1:5000 2) 1:10000